



**INNOVATIONS.
FOR YOU. BY DMG.**

Digital Workflow

Preservation

Esthetics



Useful new ideas
that make your work easier.

That's
why we're here.

What innovations does DMG have for you?

This brochure will introduce them to you, clearly and concisely. Whether you're looking for scannable precision impression material or quick and easy protection for the proximal region, we hope that our new product solutions will make your everyday work easier.

This is the aim of our over 50 years of research, development and production activities, carried out entirely in-house at our facilities in Hamburg, Germany.

Product solutions from practitioners for practitioners.

Practicing dentists are involved in the development and assessment of products and concepts at virtually every stage of our work. Initiatives such as the »Circle of Experts« and »Scientific Talk«, along with close cooperation with universities also foster exchange.

Safe and reliable.

You can rely on DMG materials »Made in Germany«. Neutral testing institutions confirm the high quality of our work for you time and again. DMG products such as Honigum, Luxatemp and LuxaCore have received top marks and continue to do so.

Honigum Pro Impressions made your way.



Stability or flowability?

The precision VPS-material makes impressions your way, on your schedule.

Its most significant benefit compared with conventional impression materials: Honigum Pro flows only when it needs to flow, i.e. under pressure. Even the smallest spaces, such as the area beneath the preparation line, are reliably reached in moist conditions. When force is not applied, Honigum Pro remains extremely stable. The material remains reliably at the application point, with no drips or run-off.

Thanks to the outstanding toughness, the impression can be safely removed without tearing. The »Snap-Set« technology ensures convenient working times.

Indications

Crown and bridge impressions, inlay and onlay impressions, all types of pickup impression (e.g., for implants), fixing of impression posts and transfer copings

Specialties

All versions of Honigum Pro are also available in a »Fast« version.

Benefits for the dentist

- Extraordinary stability AND flowability
- Optimum behavior in moist conditions
- Comfortable working time
- Extreme toughness
- Highest precision



Product video and video tutorial:

www.youtube.com/dmgdentalEN

More information at www.dmg-dental.com/honigumpro/en

Honigum Pro Scan

Impression precision made scannable.



There's nothing easier. The new scannable variants combine the superior properties of Honigum Pro with all the advantages of a digital workflow – with no need for pretreatment with scan spray or powder. You can work in your clinic as you always have: conventionally, precisely, and reliably. Without any new technology or investments, all options are available to you. In the lab, accurate impressions can be captured digitally. Possible sources of error with conventional model preparation using plaster are eliminated. Elaborate and sensitive cut-saw model preparation is not necessary. The lab's digitalized data are available to you at any time. It is also no longer necessary to store plaster models or impressions. The data permit the dentist and the lab to exchange information as early as the planning stage, regardless of location.

Take impressions easily and precisely with Honigum Pro Scan for digital processing as well.

Indications

All types of impressions such as crown and bridge impressions, inlay and onlay impressions, all types of pickup impression (e.g., for implants), fixing of impression posts and transfer copings

Specialties

Highest precision, without powder or spray

Benefits for the dentist

- Scannable: All Honigum Pro advantages for digital processing as well
- Cost effective, easy, precise – without the need to purchase expensive technology
- Reliable reproduction even of areas that are difficult to capture using an intraoral system
- Sources of error when using conventional model preparation are eliminated
- Elaborate saw-cut models are not required



O-Bite Scan

Most precise bite registration,
now also in digitized format.



O-Bite Scan combines precise bite registration with all the advantages of a digital workflow. The material is one of the hardest VPS-based impression materials. The high final hardness is particularly important for repositioning the bite registration between the models in the articulator, in order to guarantee precision-fit restorations. Softer bite registrations flex or compress, thus prohibiting a precise articulation. This can lead to unwanted bite elevations after restoration, requiring a time-consuming grinding process.

The new scannable variant means that you can also make use of the advantages of a digital workflow, without investing in equipment and without using scan powder or spray.

With O-Bite Scan you can work in your clinic as you always have: conventionally, precisely, and reliably. And then all options are available to you. The data from the accurate registration are available in digitized format at any time.

Indications

Bite registration

Specialties

Maximum break resistance and flexural strength

Benefits for the dentist

- Scannable
- Extremely high final hardness and superior break resistance
- Minimal bite resistance
- Fresh fruity orange scent
- Cost effective, easy and precise for the digital workflow – without the need to purchase expensive technology
- Easy handling, no pretreatment with powder or spray





The MiniDam protects the proximal area during treatment – quickly, easily and comfortably for the patient. The practical silicone shield ensures that the site being treated is sufficiently dry. At the same time, the gingiva is protected from the materials being used, e.g. etching gel.

The MiniDam can be placed by one person in just a few seconds and stabilizes itself – without clamps – with open access to the treatment area. That makes it extremely convenient for the dentist and very comfortable for the patient.

The material is latex-free and thus harmless for allergy sufferers. Thanks to the MiniDam, procedures in the proximal area, which also include caries infiltration for example, can thus be completed more quickly.

Indications

Silicone protection to facilitate a relatively dry environment in the proximal area between two teeth and to protect the gingiva against dental materials.

Specialties

Can be applied accurately to the treatment site and offers secure hold without clamps.

Benefits for the dentist

- Reliable local protection of the proximal region
- Quick and easy
- Comfortable for patients
- Secure hold without clamps
- Latex-free
- Assists with infiltration treatment

Available from summer 2015



Product video and video tutorial:
www.youtube.com/dmgdentalEN

More information at www.dmg-dental.com/minidam/en



This premium composite-based cementation system has been specially developed for demanding veneering procedures.

The Vitique system includes a Cementation Kit with cements and try-in pastes, as well as an Accessory Kit with etching gel, silane and bonding agent. The key benefits of the materials include flexible shade selection, an exact match between try-in and cement, maximum color stability and easy excess removal. The Vitique system also provides you with comprehensive support with practice-related services. The product comes with a comprehensive illustrated user guide that clearly presents the entire treatment procedure. Moreover, a special expert hotline has been set up through which renowned veneering specialists offer practical assistance.*

Indications

Definitive, highly esthetic cementation of ceramic or composite veneers

Specialties

Complete system for esthetic cementation of veneers comprised of a Cementation Kit and Accessory Kit

Benefits for the dentist

- Complete, superior-quality system for veneer cementation
- Outstanding shade reproduction and shade selection
- Exact shade match between try-in and cement
- Particularly easy removal of excess material
- High shade stability
- Comprehensively illustrated user guide included
- Additional application support and value-added services: Video tutorials, photo course, expert hotline...



Product video and video tutorial:
www.youtube.com/dmgdentalEN

All steps at a glance: Comprehensive illustrated step-by-step user guide with 6 step-by-step sheets included with the packaging.

More information at www.dmg-dental.com/vitique/en

*Please note that our experts cannot provide diagnostic and therapeutic advice, but are only available for advising on the optimal esthetic use of veneers.



This flowable composite is not only self-adhesive, but also self-etching, making the practice workflow considerably easier. You are spared the etching and bonding steps and the associated time expenditure. Potential sources of error are minimized as well, benefits that conventionally adhesively-cemented flowables do not offer.

Constic offers excellent adhesion to enamel and dentine and helps prevent post-operative sensitivity.

The flowable composite can be used for small class I and III restorations, as well as for base lining and for fissure sealings. Class V restorations and the blocking out of undercuts are also quick and easy. The material is simply rubbed in using a brush and light-cured.

Indications

Small class I and III restorations, class V restorations, (extended) sealing of fissures, base lining of class I-II restorations, blocking out and filling of undercuts in cavities, small occlusal cavities in deciduous teeth

Specialties

Very good adhesion to enamel and dentine without a separate etching and bonding step

Benefits for the dentist

- Self-etching and self-adhesive
- Faster, easier and more efficient treatment procedure
- Radiopacity and fluorescence similar to that of teeth
- Prevents postoperative sensitivities
- Minimizes potential sources of mistakes

Product video and video tutorial: www.youtube.com/constic

All steps at a glance:
Illustrated step-by-step sheet included with the packaging.

More information at www.dmg-dental.com/constic/en



PACKAGING

Honigum Pro

Package:	4 Cartridges @ 25 ml, 10 Automix-Tips, 10 Intraoral-Tips	
	Honigum Pro-Light	REF 929539
	Honigum Pro-Light Fast	REF 929639
	2 Cartridges @ 50 ml, 10 Automix-Tips, 10 Intraoral-Tips	
	Honigum Pro-Light	REF 929831
	Honigum Pro-Light Fast	REF 929832
	2 Cartridges @ 50 ml, 4 Automix-Tips, 4 Intraoral-Tips	
	Honigum Pro-Mono	REF 929830
	Honigum Pro-Mono Fast	REF 929976
	2 Cartridges @ 50 ml, 4 Automix-Tips	
	Honigum Pro-Heavy	REF 929766
	Honigum Pro-Heavy Fast	REF 929838
	1 Cartridge @ 380 ml, 10 MixStar-Tips	
	Honigum Pro-Mono	REF 929568
	Honigum Pro-Mono Fast	REF 929978
	Honigum Pro-Heavy	REF 929537
	Honigum Pro-Heavy Fast	REF 929638
	Honigum Pro-Putty	REF 929374

Honigum Pro-Heavy Scan

Package:	2 Cartridges @ 50 ml, 4 Automix-Tips	REF 929983
	1 Cartridge @ 380 ml, 10 MixStar-Tips	REF 929981

Honigum Pro-Light Scan

Package:	2 Cartridges @ 50 ml, 10 Automix-Tips	REF 929984
----------	---------------------------------------	------------

O-Bite Scan

Package:	2 Cartridges @ 50 ml, 12 Automix-Tips	REF 919017
----------	---------------------------------------	------------

Constic

Package:	2 Syringes @ 2 g composite, 1 Brush holder, 20 Clip-on brushes, 20 Luer-Lock-Tips Silver			
	Shade A1	REF 220700	Shade A3.5	REF 220703
	Shade A2	REF 220701	Shade B1	REF 220704
	Shade A3	REF 220702	Shade Opak-White	REF 220705

Vitique

Cementation System

	3 Syringes @ 6 g ESTHETIC RESIN CEMENT	
	3 Syringes @ 3.9 g TRY-IN-PASTE in the shades A2.5; B1; TR	
	1 Syringe @ 2 ml Etching Gel	
	1 Bottle @ 5 ml TECO	
	1 Bottle @ 3 ml Vitique Silane Adhesive	
	1 Bottle @ 3 ml Vitique Silane Activator	
	Accessories	REF 213080

MiniDam

Package:	20 Pieces	REF 220381
----------	-----------	------------

DMG

Chemisch-Pharmazeutische Fabrik GmbH
Elbgaustraße 248 22547 Hamburg Germany
Fon: +49. (0) 40. 84 006-0 Fax: +49. (0) 40. 84 006-222
info@dmg-dental.com www.dmg-dental.com
www.facebook.com/dmgdental