

Pro AF Baby Gold




Indications


- 1. Pediatric Rotary Endodontics
- Adult Rotary endodontics in condition where there is limited accessibility eg. Third molar root canal treatment, OSMF

Salient Features 🛛 🚤

- Short 17mm file improving safety with comfort to both dentist and the patient
- Advanced NiTi CM wire for better canal centricity unlike NiTi file.
- Heat Treatment resulting nearly no separation
- A versatile rotary file system suitable for conservative preparation of all types of canals.
- Improved Shaping of canals with sequential combination of 4% and 6% taper files



Accurate tip grinding Non cut tip




Sharp Edge, Super high cutting efficiency



Well Polished surface by advanced patented technology

Instrumentation Protocol:

Pro AF Baby Gold



Posterior


Anterior






Free Autoclavable Endo Box

Cross Section of Root under CBCT



Unprepared Flattened canal



Conservative Preparation (4% taper file)

Aggressive

Preparation

(4% taper file) (6% taper file) Unmatched cleaning efficiency with conservative canal preparation in ribbon shaped canals as compared to variable taper rotary file

Cases -by Dr. Shaymaa Moustafa (Egypt)

